

Water Flosser Tips



**Classic
Jet Tip**

For general use.
Cleans deep
between teeth and
below the gumline.

Replace every
6 months



**Pik Pocket®
Tip**

For periodontal
pockets and
furcations.
A soft rubber tip
designed for gentle,
low-pressure delivery
of therapeutic rinses
below the gumline.
Clinically proven
to access 90%
of a 6mm pocket.

Replace every
3 months



**Orthodontic
Tip**

For orthodontic
appliances and
dental braces.
The innovative
tapered brush
removes hard-to-
reach plaque around
orthodontic brackets
and wires.

Replace every
3 months



**Plaque Seeker®
Tip**

For implants, crowns,
bridges and other
restorations.
Three thin tufts of
bristles help to access
stubborn plaque.

Replace every
3 months



**Tongue
Cleaner Tip**

For fresh breath.
With a gentle
scraping and flushing
action, the Waterpik™
Tongue Cleaner
removes bacteria and
sulfur compounds
that reside on the
tongue and cause
bad breath.

Replace every
6 months



**Toothbrush
Tip**

Stationary toothbrush
tip allows water
to pulsate through
the brush to clean
and rinse the tooth
surface. Designed
to let you brush and
water floss at the
same time.
Use with or without
toothpaste.

Replace every
3 months

waterpik™
waterflosser™

Join our Facebook Community:



@waterpikprofessionalau

www.waterpikshop.com.au



RECOMMENDED BRAND
BY AUSTRALIAN DENTISTS*

*Based on a survey of 78 Australian dentists conducted
by a leading research company in February 2022.